

## Introduction

**Product Name:** Infant SpO<sub>2</sub> Sensor with Disposable wrap

**Model Number:** 15040072

**Device Compatibility:** CMI PC-66 series ONLY

**Components:** It consists of 3 parts.

1. Y-style SpO<sub>2</sub> Sensor
2. Disposable Foot Wrap
3. Ankle Wrap

**Recommended Patients:** Infants and toddlers under 25 lbs.

## Instructions for use




Figure A

1. For proper placement on either foot, place the sensors on the outside of the foot behind the pinky toe. Make sure the sensor touch the skin closely, then secure the foot wrap with Velcro (see Figure A and B). Do not over-tighten.




Figure B

Note: Improper sensor placement may result in difficulty acquiring signals or inaccurate results

2. Use the ankle wrap to secure the sensor cable on the ankle or leg (see Figure C). Do not over-tighten.




Figure C Right foot view

Note: The device display may show inaccurate reading or “Check Probe” due to motion artifact induced by the action of sensor placement. Power off and re-start the device after sensor placement to obtain accurate readings.

## Change Disposable Wrap




Figure D Y-style sensor



Figure E Foot wrap

Press the light source and receiver into the slots of the foot wrap respectively. Avoid pulling the silicon part of the sensor with force. Doing so may damage the sensor integrity and cause sensor malfunction. Please note that the light source, receiver and the Velcro should be facing the same side of the wrap (see Figure F).




Figure F

## **Intended Use**

This sensor is intended to be used with a compatible CMI Health brand pulse oximeter for measuring the functional oxygen saturation (SpO<sub>2</sub>) and pulse rate of a specified patient type.

## **Attentions**

- 🔔 The operation of this sensor can only be performed by trained personnel.
- 🔔 ALWAYS wrap and secure the sensor before switching on the CMI Health Oximeter device. (Wrapping with the device on may give inaccurate readings)
- 🔔 If the sensor is wrapped too tightly, circulation may be blocked, which leads to discomfort and inaccurate readings.
- 🔔 Make sure the light source and receiver of the Y-style sensor are placed near the edge of the instep, and wrapping is done with proper tightness.
- 🔔 If the sensor does not provide reliable pulse signal, it may be incorrectly positioned. If such situation occurs, reposition the sensor on the foot until a reliable pulse signal can be detected.
- 🔔 Strong surrounding light sources, such as fluorescent light, ruby lamp, infrared heating lamp, and direct sunlight, may cause inaccurate readings.
- 🔔 Excessive patient movement and the extremely strong electromagnetic interference may cause unstable signals and inaccurate readings.

## **Warnings**

- ⚠️ Do not alter or modify the sensor. Alterations and modifications may affect performance or accuracy.
- ⚠️ This sensor should be used together with the compatible oximetry device, otherwise the sensor may not work or the reading may be inaccurate.
- ⚠️ Although the biocompatibility evaluation has been performed on this sensor, some exceptional allergic patients may still cause anaphylaxis. Do not apply this sensor to those who has anaphylaxis.

- ⚠️ Change the measuring site every 2 or 3 hours. When the ambient temperature is over 35°C, change the measuring site every 2 hours. When the ambient temperature is over 37°C, STOP using this sensor immediately since long time measurement may cause serious scalding or burn injury.
- ⚠️ The measuring site must be examined more carefully for patients with special conditions. Do not place the sensor on the site with edema or fragile tissue.
- ⚠️ Misapplication of the sensor with excessive pressure for prolonged periods can induce pressure injury.
- ⚠️ Check the integrity of the sensor before use, discard and replace the sensor if it is damaged.

Note: This sensor is compatible only with devices sold by CMI Health or the CMI Store.

Other information about this sensor, please refer to the user manual of its compatible device.

## **Specifications**

SpO<sub>2</sub> measuring range: 35%~100%

SpO<sub>2</sub> measuring accuracy: A<sub>rms</sub> value (defined in ISO 9919 /ISO 80601-2-61) is not greater than 3% in the range of 70%~100%.

Pulse Rate measuring range: 30bpm~250bpm

Pulse Rate measuring accuracy: ±2bpm or ±2%,

whichever is greater.

Wavelength: Red light: 663nm, Infrared light: 890nm



Do not litter at will



Refer to the accompanying documents

---

**CMI Health, Inc.**

**5975 Shiloh Road, Suite 114**

**Alpharetta, GA 30005, USA**

**Toll-free: 888-985-1125**

**E-mail: [info@cmihealth.com](mailto:info@cmihealth.com)**

**Website: [www.cmihealth.com](http://www.cmihealth.com)**